



H i g h T e c h n o l o g y > Q u a l i t y > P r o f e s s i o n a l i t y > E x p e r i e n c e

EPI-C plus®



TECHNICAL DATA:

Maximum Energy	760J
Fluence	From 18J/cm ² to 46J/cm ²
OPE® Filters/Lamps	390 nm Phototherapy
	550 nm Skin Rejuvenation
	570 nm Hair Removal (light skins)
	630 nm Hair Removal (dark skins)
	S.M. Stretch Marks
	LIPS Lips
YOUNG AGAIN®	RED: Stimulates production of new collagen
	YELLOW: Specific action on lymphatic system
	BLUE: Skin purification
Pulse Length	From 10 to 30 milliseconds
Lamps Spot Area	2,5 cm x 5 cm = 12,5 cm ²
Control Panel	PC Touch Screen
Supply Voltage	230 Vac 50 Hz - 4A Other voltage on request (e.g. 110 Vac 60 Hz 8A)
Cooling System (Patented)	Forced air circulation inside the handpiece
Weight	Kg 25
Dimensions	30 cm x 47cm x 45 cm



EPI-C^{plus}®

The synergy of Polychromatic and Monochromatic Light

The Technologies:

OPE® Polychromatic High Intensity Light

YOUNG AGAIN® Monochromatic Biostimulating Light

The Treatments:

FACE

- Antiaging
- Impure skin
- Greasy skin
- Filling out lips
- Photo-Aging
- Chrono-Aging
- Toning Face/Neck/Décolleté
- Permanent Epilation

HEAD

- Hair strengthening

BODY

- Tissue oxygenation
- Breast
- Buttocks
- Stretch marks
- Permanent epilation
- Hands rejuvenation